3328R  
Lot 011~

**MBL**  
A JSR Life Sciences  
Company

Printed March 5, 2024  
Version 7.0

For Research Use Only. Not for use in diagnostic procedures.

## **DDDDK-tagged Protein PURIFICATION GEL**

(MoAb. clone FLA-1GS)

CODE No. 3328R

*PURIFICATION to Maintain Protein Activity  
from eukaryote cell lysate and culture supernatant*

## **Product Description**

The ability to isolate and study a purified protein lies at the heart of modern biochemistry. Researchers in many fields require highly purified, active proteins for studies involving signaling pathways, enzymology, receptor binding, DNA binding, post-transcriptional modifications, and much more. Thus, choosing a method of purification is an important aspect in maintaining protein structure and function.

Recombinant tagged protein purification methods and kits are now widely recognized. The DDDDK epitope tag (DYKDDDDK) has been widely used as a multi-purpose tag, and anti-DDDDK-tag antibodies are optimally suited for identifying, detecting, purifying, and monitoring the expression levels of recombinant DDDDK-tagged proteins.

MBL's DDDDK-tagged Protein PURIFICATION GEL is designed for the isolation of DDDDK-tagged protein from cell culture supernatants and cell lysate under neutral pH condition. Severe conditions such as acidic or alkaline elution denature protein structure. However, a neutral pH elution can preserve protein activity and native conformation. Therefore, MBL has developed the Anti-DDDDK-tag Gel to purify DDDDK-tagged proteins quickly and efficiently at neutral pH to maintain the activity and conformation of purified proteins. The elution of DDDDK-tagged proteins from the Gel is achieved by the addition of the DDDDK-tag peptide (DYKDDDDK). As the DDDDK-tag peptide competes with DDDDK-tagged proteins on the Gel, the purified proteins do not lose the protein activity.

## **Components**

	Quantity
	CODE No. 3328R
<b><u>Anti-DDDDK-tag Gel</u></b>	<u>5 mL × 1 vial</u>

50% slurry: 5 mL Gel in 10 mL total volume in PBS with 0.1% ProClin 150 as preservative.

Elution Peptide (DYKDDDDK) is NOT included.

It is available as CODE No. 3325-205 (1 mg x 5 vials).

## **Product Capacity**

The purification capacity of the Anti-DDDDK-tag Gel varies depending upon the DDDDK-tagged protein. For example, 1 mL of anti-DDDDK tag gel purified 0.3 mg of DDDDK-tagged protein.

## **Storage**

Store for up to 1 year from date of receipt at 2-8°C. Do not freeze.

## Material Preparation

Prepare the following reagents before affinity purification.

Lysis buffer : Suitable Lysis buffer varies with a kind of the DDDDK-tagged protein  
(see **Additional Information**).

10-50 mM Tris-HCl (pH 7.5)  
100-300 mM NaCl  
1-5 mM EDTA  
1% NP-40 or Triton X-100  
if necessary add Protease Inhibitor Cocktail  
(e.g. SIGMA: code P8340, PIERCE: code 78415)

Washing buffer : Suitable Washing buffer varies with a kind of the DDDDK-tagged protein  
(see **Additional Information**).

10-50 mM Tris-HCl (pH 7.5) or HEPES-KOH (pH 7.5)  
300-500 mM NaCl  
0.1% NP-40 or Triton X-100


To prevent non-specific protein binding to the gel, salt concentration should be 300-500 mM NaCl.

Elution buffer : 0.1 mg/mL DDDDK-tag peptide in Washing buffer.

Regeneration buffer : 0.17 M Glycine-HCl, pH 2.3

Column storage buffer : PBS/Preservative (e.g. 0.1% ProClin 150 and 0.09% NaN<sub>3</sub>)

## Procedure Summary

The descriptions of the following protocols are examples. Each user should determine the appropriate condition.

## **Protocols**

- Notes:
1. The Gel is optimized for the native conditions only. It is not recommended to use under denaturing conditions, nor for purification of aggregated, unstable, and insoluble protein (e.g. inclusion bodies). Proteins solubilized with such as 6 M Guanidine-HCl or 8 M Urea cannot be purified using this Gel (see **Additional Information**).
  2. Cellular debris and particulate matter must be removed prior to purification. The protein extract should be centrifuged ( $10,000-20,000 \times g$  for 15 minutes) and filtered with a  $0.45 \mu\text{m}$  filter to remove any remaining cells and particulates.
  3. Highly viscous samples containing chromosomal DNA or RNA should be sonicated or treated with nuclease to reduce viscosity.

### A. Column preparation

1. Place the empty chromatography column (e.g., PIERCE, code 29920) on rack or stand.
2. Rinse the column with Washing buffer.
3. Resuspend Anti-DDDDK-tag Gel by tapping and inverting the vial several times immediately before dispensing. Don't vortex.
4. Transfer the desired volume to the column. Allow the column to drain.  
Do not allow the column to dry out.
5. Wash the column with 10 bed volumes of Washing buffer.

### B. Column binding

1. Apply the lysate on the top of the column under gravity flow.

Note: The binding efficiency to Anti-DDDDK-tag Gel depends on the conditions such as a kind of the DDDDK-tagged protein, sample flow rate, and temperature. If the coupling efficiency has been low, put a sample through the column several times, or incubate the lysate with the Gel in batch.

2. Collect the flow-through into clean collection tubes.
3. Wash the column with Washing buffer until the OD<sub>280</sub> is  $<0.01$ .

### C. Elution of the DDDDK-tagged protein

1. Prepare 8 bed volumes of Elution buffer.
2. Remove the bottom cap, and pour off the excess liquid.

3. To exchange the buffer in the column, add 1 bed volume of Elution buffer, and drain the 1 bed volume of the buffer.
4. Cap the bottom, and add 1 bed volume of Elution buffer so that the Gel does not dry.
5. Allow the column to stand at 4°C for 5 minutes. (As the DDDDK-tag peptide competes with DDDDK-tagged protein on the Gel by this incubation, DDDDK-tagged protein is eluted from the Gel.)
6. Remove the bottom cap, and collect the eluate into clean collection tubes (Fraction 1).
7. Add 1 bed volume of Elution buffer, collect the eluate into clean collection tubes (Fraction 2).
8. Repeat the elution and collection step (Step 6 and 7) another 5 times. (Fraction 1-7 contains DDDDK-tagged protein.)

#### D. Regeneration and storage

1. Wash the column with 10 bed volumes of Regeneration buffer.
2. Immediately wash the column with 10 bed volumes of Column storage buffer.
3. Store at 2-8°C in 2 bed volumes of Column storage buffer.


Note: Poured columns containing Anti-DDDDK-tag Gel may be used at least 10 times, depending on the usage conditions.

#### **Related Products:**

For more information, please visit our website at <https://ruo.mbl.co.jp/>.

## Example of Purification Results 1

### Purification of N-terminus DDDDK-tagged $\beta$ -galactosidase


#### SDS-PAGE (Coomassie Brilliant Blue Staining)

Human embryonic kidney cells (293T) were transfected with pcDNA-DDDDK-tagged  $\beta$ -galactosidase and cultured for 60 hours. Cells were then lysed in the Lysis buffer (10 mL/100-mm dish x 5) and purified according to the preceding protocol. Column bed volume was 0.25 mL. Elution was carried out with 2 mL of 0.1 mg/mL DDDDK-tag peptide. Each fraction was 0.25 mL.

## Example of Purification Results 2

### Purification of C-terminus DDDDK-tagged GFP from *E. coli* lysate

#### SDS-PAGE (Coomassie Brilliant Blue Staining)

DDDDK-tagged GFP plasmid was transformed into *E. coli* BL21(DE3)RIL and cultured in 5 mL LB media at 27°C for 12 hours. Cells were then lysed in the 5 mL Lysis buffer and purified according to the preceding protocol. Column bed volume was 1 mL. Elution was carried out with 8 mL of 0.1 mg/mL DDDDK-tag peptide. Each fraction was 1 mL.

## Additional Information

Several reagents were examined whether they were suitable for use with the DDDDK-tagged Protein PURIFICATION GEL. For example, RIPA buffer could be used for preparation of cell lysate. The results are listed below.

### Chaotropic agents

<b>Urea</b>	<b>1 M</b>	<b>Yes</b>
<b>Guanidine-HCl</b>	<b>1 M</b>	<b>No</b>

### Reducing agents

<b>DTT</b>	<b>10 mM</b>	<b>Yes</b>
<b>2-Mercaptoethanol</b>	<b>10 mM</b>	<b>Yes</b>

### Surfactants

<b>Nonionic</b>	<b>Tween-20</b>	<b>5%</b>	<b>Yes</b>
	<b>TritonX-100</b>	<b>5%</b>	<b>Yes</b>
	<b>NP40</b>	<b>1%</b>	<b>Yes</b>
	<b>Digitonin</b>	<b>1%</b>	<b>Yes</b>
	<b>n-Octyl-<math>\beta</math>-D-gulcoside</b>	<b>1%</b>	<b>Yes</b>
<b>Zwitterionic</b>	<b>CHAPS</b>	<b>1%</b>	<b>Yes</b>
	<b>CHAPSO</b>	<b>1%</b>	<b>Yes</b>
<b>Anionic</b>	<b>SDS</b>	<b>0.1%</b>	<b>Yes</b>
	<b>Sodium Deoxycholate</b>	<b>0.5%</b>	<b>Yes</b>

### Others

<b>NaCl</b>	<b>1 M</b>	<b>Yes</b>
<b>Glycerol</b>	<b>10%</b>	<b>Yes</b>
<b>EDTA</b>	<b>10 mM</b>	<b>Yes</b>

The “Yes” indicates the reagents can be used in the Lysis buffer for this Gel up to the indicated concentration. The “No” indicates the reagents cannot be used in the Lysis buffer for this Gel at the indicated concentration.

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