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Code No.
8469

Quantity
600 tests

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Quantity
600 tests

PRODUCT DESCRIPTION

Histostar™ DAB Substrate is designed for use on immunohistochemical staining of tissue sections. It produces a brown colored deposit at the site of antigen upon reaction with peroxidase.

KIT COMPONENTS

Chromogenic substrate (Solution A)	3,3'-Diaminobenzidine (DAB) • 4HCl	3 mL x 1 bottle
Substrate buffer (Solution B)		3 mL x 1 bottle
Chromogenic reagent (Solution C)	0.6% H ₂ O ₂	3 mL x 1 bottle

STORAGE and STABILITY

Store at 2–8°C. Stable for minimum of 18 months from manufacturing date when stored at 2–8°C. The expiration date is printed on the reagent vial.

The descriptions of the following protocols are examples. Each user should determine the appropriate condition.

PROTOCOL

1) Add 1 drop (or 25 µL) of the solution A and solution B into 1 mL of distilled water. Mix well, then add 1 drop (or 25 µL) of solution C, and mix again. Protect this prepared solution from light during storage, and apply the solution within 30 minutes.

*11 slides are stainable, when apply the solution 100 µL each.

2) Incubate the specimen, which has been reacted with peroxidase conjugate, with 1–2 drops of Histostar™ DAB Substrate Solution at room temperature (15–25°C) for 5–20 minutes. Monitor color development with a microscope. The incubation time will depend on the conditions. After adequate coloration, wash the slide with tap water in order to stop color development.

PRECAUTIONS

- 1) 3,3'-Diaminobenzidine (DAB) • 4HCl have carcinogenicity therefore must be handled with care. Always wear gloves.
- 2) Chromogenic reagent (Solution C) containing 0.6% H₂O₂ therefore, it must be handled with care.

PRECAUTIONS OF DISPOSE

This substrate containing 3,3'-Diaminobenzidine (DAB), therefore, it must be disposed following a proper method outlined by each lab's disposal protocol.

INTENDED USE

For Research Use Only. Not for use in diagnostic procedures.

RELATED PRODUCTS

Please visit our web site <https://ruo.mbl.co.jp/>