



Non-Radioisotopic Kit for Measuring CK2 Activity

CycLex CK2 Kinase Assay/Inhibitor Screening Kit

Cat# CY-1170

Intended Use.....	1
Storage.....	1
Introduction.....	2
Principle of the Assay.....	3
Materials Provided.....	4
Materials Required but not Provided.....	4
Precautions and Recommendations.....	5
Detailed Protocol.....	6-9
Evaluation of Results.....	10
Assay Characteristics.....	10
Troubleshooting.....	10
Reagent Stability.....	10
Example of Test Results.....	11-12
References.....	13

Intended Use

The MBL Research Product **CycLex CK2 Kinase Assay/Inhibitor Screening Kit** designed to measure the activities of purified CK2 for the rapid and sensitive evaluation of inhibitors or activators. The phospho-serine specific monoclonal antibody used in this assay kit has been demonstrated to recognize the phospho-serine46 residue in p53, which is phosphorylated by CK2 in vitro.

Applications of this kit include:

- 1) Screening inhibitors or activators of CK2.
- 2) Detecting the effects of pharmacological agents on CK2 activity.

This assay kit is for research use only and not for use in diagnostic or therapeutic procedures.

Storage

- Upon receipt store all components at 4°C.
- Don't expose reagents to excessive light.



Introduction

Protein kinase Casein kinase II (CK2) is a ubiquitous and pleiotropic seryl/threonyl protein kinase, which appears to interact with different signaling pathways and therefore represents the prototype of a multifunctional protein kinase. The holoenzyme is generally composed of two catalytic (alpha and/or alpha') and two regulatory (beta) subunits⁽¹⁻³⁾, but the free alpha/alpha' subunits are catalytically active by themselves. Although the beta subunits deeply affect many properties of CK2, both the isolated catalytic subunits and the holoenzyme are constitutively active. The enzyme is highly expressed in most cancers⁽⁴⁾ and this higher expression has been tentatively correlated with the involvement of CK2 in the promotion of specific phases of the cell cycle⁽⁵⁾. Unlike the majority of protein kinases, which are tightly regulated enzymes, CK2 is endowed with high constitutive activity, a feature that is suspected to underlie its oncogenic potential^(6, 7) and possible implication in viral infections. This makes CK2 an attractive target for anti-neoplastic and antiviral drugs.

Experimental studies suggest that dysregulated expression of the alpha subunit of CK2 imparts an oncogenic potential in the cells such that in cooperation with certain oncogenes^(8, 9), it produces a profound enhancement of the tumor phenotype. Recent studies have provided evidence that overexpression of CK2 in tumor cells is not simply a reflection of tumor cell proliferation alone but additionally may reflect the pathobiological characteristics of the tumor. Of considerable interest is the possibility that CK2 dysregulation in tumors may influence the apoptotic activity in those cells⁽¹⁰⁻¹²⁾. Approaches to interfering with the CK2 signal may provide a useful means for inducing tumor cell death⁽¹³⁾.

Measurement of CK2 activity

The protocol generally regarded as most sensitive for the quantitative measurement of CK2 activity involves incubation of the CK2 sample with substrate, either a natural or synthetic polypeptide (such as CK2 substrate peptide; RRRDDDSDDD), in the presence of Mg^{2+} and ^{32}P -labeled ATP. The reaction is terminated by "spotting" a sample onto a phosphocellulose P81 filter paper disc, followed by washing extensively to remove unincorporated radiolabel and the incorporated radioactivity on P81 filter is counted. While sensitive, this method is labor-intensive, generates hazardous radioactive waste, and depends on a radioisotope of short half-life. It is particularly unsuitable when kinase assays are only performed on an infrequent basis. The MBL Research Product **CycLex CK2 Kinase Assay/Inhibitor Screening Kit** uses a peroxidase coupled anti-phospho-p53 serine46 monoclonal antibody as a reporter molecule in a 96-well ELISA format. This assay provides a non-isotopic, sensitive and specific method to detect CK2 activity.




Principle of the Assay

The MBL Research Product **CycLex CK2 Kinase Assay/Inhibitor Screening Kit** is a single-site, semi-quantitative immunoassay for CK2 activity. Plates are pre-coated with a substrate corresponding to recombinant p53, which contains a serine residue that are phosphorylated by CK2 (Casein kinase II).

The detector antibody specifically detects only the phosphorylated form of serine46 on p53. The MBL Research Product **CycLex CK2 Kinase Assay/Inhibitor Screening Kit** can be used to study the kinetics of a purified or partially purified CK2 as well as to screening these kinases inhibitor. To perform the test, the sample is diluted in Kinase Buffer, pipetted into the wells and allowed to phosphorylate the bound substrate in the presence of Mg^{2+} and ATP. The amount of phosphorylated substrate is measured by binding it with a horseradish peroxidase conjugate of TK-4D4, an anti-phospho-p53 serine46 specific antibody, which then catalyzes the conversion of the chromogenic substrate tetra-methylbenzidine (TMB) from a colorless solution to a blue solution (or yellow after the addition of stopping reagent). The color is quantified by spectrophotometry and reflects the relative amount of CK2 activity in the sample. For kinetic analysis, the sample containing CK2 is added to the wells in a similar fashion and at varying times the reaction is stopped by the addition of a chelator, sodium ethylenediaminetetraacetate (EDTA) and the amount of phosphorylated substrate determined as before.

The MBL Research Product **CycLex CK2 Kinase Assay/Inhibitor Screening Kit** is designed to accurately determine the presence and relative amount of CK2 activity in purification column fractions, and to determine non-isotopic kinetic analysis of CK2 activity. Careful attention to methods of chromatography and the assay protocol will provide the investigator with a reliable tool for the evaluation of CK2 activity.

Summary of Procedure





Materials Provided

All samples and standards should be assayed in duplicate. The following components are supplied and are sufficient for the one 96-well microtiter plate kit.

Microplate: One microplate supplied ready to use, with 96 wells (12 strips of 8-wells) in a foil, zip-lock bag with a desiccant pack. Wells are coated with recombinant p53 N-terminus (1-99a.a.) as substrate of CK2.

10X Wash Buffer: One bottle containing 100 mL of 10X buffer containing Tween[®]-20

Kinase Buffer: One bottle containing 20 mL of 1X buffer; used for Kinase Reaction Buffer and sample dilution.

20X ATP: One vial of lyophilized ATP Na₂ salt.

HRP conjugated Detection Antibody: One vial containing 12 mL of HRP (horseradish peroxidase) conjugated anti-phospho-p53 S46 monoclonal antibody (TK-4D4). Ready to use.

Substrate Reagent: One bottle containing 20 mL of the chromogenic substrate, tetra-methylbenzidine (TMB). Ready to use.

Stop Solution: One bottle containing 20 mL of 1 N H₂SO₄. Ready to use.

Materials Required but not Provided

- **CK2 (alpha/beta) Positive Controls:** Available from MBL (CK2 (alpha/beta): Cat# CY-E1170-1).
- **(Optional) 10X Heparin (1 µg/mL):** As a CK2 inhibitor, available from Sigma, Cat# H-4784. 100 µg/mL stock solution (H₂O) diluted 1:100 in Kinase Buffer.
- **Pipettors:** 2-20 µL, 20-200 µL and 200-1,000 µL precision pipettors with disposable tips.
- **Wash bottle or multichannel dispenser** for plate washing.
- **Microcentrifuge and tubes** for sample preparation.
- **Vortex mixer**
- **(Optional) Microplate washer:** Manual washing is possible but not preferable.
- **Plate reader:** capable of measuring absorbance in 96-well plates at dual wavelengths of 450 nm/540 nm. Dual wavelengths of 450/550 or 450/595 nm can also be used. The plate can also be read at a single wavelength of 450 nm, which will give a somewhat higher reading.
- **(Optional) Software package facilitating data generation and analysis**
- **500 or 1,000 mL graduated cylinder**
- **Reagent reservoirs**
- **Deionized water of the highest quality**
- **Disposable paper towels**



Precautions and Recommendations

- **Although we suggest to conduct experiments as outlined below, the optimal experimental conditions will vary depending on the parameters being investigated, and must be determined by the individual user.**
- Allow all the components to come to room temperature before use.
- All microplate strips that are not immediately required should be returned to the zip-lock pouch, which must be carefully resealed to avoid moisture absorption.
- Do not use kit components beyond the indicated kit expiration date.
- Use only the microtiter wells provided with the kit.
- Rinse all detergent residue from glassware.
- Use deionized water of the highest quality.
- Do not mix reagents from different kits.
- The buffers and reagents in this kit may contain preservatives or other chemicals. Care should be taken to avoid direct contact with these reagents.
- Do not mouth pipette or ingest any of the reagents.
- Do not smoke, eat, or drink when performing the assay or in areas where samples or reagents are handled.
- Dispose of tetra-methylbenzidine (TMB) containing solutions in compliance with local regulations.
- Avoid contact with the acidic Stop Solution and Substrate Solution, which contains hydrogen peroxide.
- Wear gloves and eye protection when handling immunodiagnostic materials and samples of human origin, and these reagents. In case of contact with the Stop Solution and the Substrate Solution, wash skin thoroughly with water and seek medical attention, when necessary.
- **Biological samples may be contaminated with infectious agents. Do not ingest, expose to open wounds or breathe aerosols. Wear protective gloves and dispose of biological samples properly.**
- **CAUTION: Sulfuric Acid is a strong acid. Wear disposable gloves and eye protection when handling Stop Solution.**



Detailed Protocol

The MBL Research Product **CycLex CK2 Kinase Assay/Inhibitor Screening Kit** is provided with removable strips of wells so the assay can be carried out on separate occasions using only the number of strips required for the particular determination. Since experimental conditions may vary, an aliquot of the CK2 Positive Controls (See the section “Materials Required but not Provided” above) should be included in each assay. Disposable pipette tips and reagent troughs should be used for all liquid transfers to avoid cross-contamination of reagents or samples.

Preparation of Working Solutions

1. Prepare a working solution of **Wash Buffer** by adding 100 mL of the **10X Wash Buffer** (provided) to 900 mL of deionized (distilled) water (ddH₂O). Mix well. Store at 4°C for two weeks or -20°C for long-term storage.
2. Prepare **20X ATP Solution** by adding **0.8 mL** of ddH₂O to the vial of **20X ATP** (provided, lyophilized). Mix gently until dissolved. The final concentration of the **20X ATP Solution** should be **2.5 mM**. Store the solution in small aliquots (e.g. 100 µL) at -20°C.
3. Prepare **Kinase Reaction Buffer** by mixing following reagents.

	96 assays	10 assays	1 assay
Kinase Buffer (provided)	9.5 mL	950 µL	95 µL
20X ATP Solution	0.5 mL	50 µL	5 µL
Total	10 mL	1000 µL	100 µL

You will need 80-90 µL of Kinase Reaction Buffer per assay well. Mix well. Discard any unused Kinase Reaction Buffer after use.

Standard Assay

1. Remove the appropriate number of microtiter wells from the foil pouch and place them into the well holder. Return any unused wells to the foil pouch, re-fold, seal with tape and store at 4°C.
2. Prepare all samples (diluted with Kinase Buffer as needed). All samples should be assayed in duplicate.
3. Add **10 µL** of diluted enzyme samples to the wells of the assay plate on ice. Duplicate wells containing 10 µL of CK2 Positive Controls (2 m units/µL) should be included in each assay as a positive control for phosphorylation.
4. Begin the kinase reaction by addition of **90 µL** of **Kinase Reaction buffer** per well, cover with plate sealer, and incubate at 30°C for 30 minutes.
5. Wash wells five times with Wash Buffer making sure each well is filled completely. Remove residual Wash Buffer by gentle tapping or aspiration.



For Research Use Only, Not for use in diagnostic procedures

6. Pipette **100 μ L** of **HRP conjugated Detection Antibody** into each well, cover with a plate sealer and incubate **at room temperature (ca.25°C) for 30 minutes**. Discard any unused conjugate.
7. Wash wells five times as same as in step 5.
8. Add **100 μ L** of **Substrate Reagent** to each well and incubate **at room temperature (ca.25°C) for 5–15 minutes**.
9. Add **100 μ L** of **Stop Solution** to each well in the same order as the previously added Substrate Reagent.
10. Measure absorbance in each well using a spectrophotometric plate reader at dual wavelengths of 450/540 nm. Dual wavelengths of 450/550 or 450/595 nm can also be used. Read the plate at 450 nm if only a single wavelength can be used. Wells must be read within 30 minutes of adding the Stop Solution.

Kinetic Assay

1. Remove the appropriate number of microtiter wells from the foil pouch and place them into the well holder. Return any unused wells to the foil pouch, refold, seal with tape and store at 4°C.
2. Prepare all samples (diluted with Kinase Buffer as needed). All samples should be assayed in duplicate.
3. Add **10 μ L** of diluted enzyme samples to the wells of the assay plate on ice. Duplicate wells containing 10 μ L of CK2 Positive Controls (2 m units/ μ L) should be included in each assay as a positive control for phosphorylation.
4. Begin kinase reaction by addition of **90 μ L** of **Kinase Reaction Buffer** in duplicate per well in timed intervals (suggested interval is 5 minutes but should be individually determined for each system). After the final addition, incubate **at 30°C for 20 minutes**.
5. Stop the reaction by flicking out the contents. (Alternatively, the reaction may be terminated by the addition of 150 μ L 0.1 M Na EDTA, pH 8.0 to each well).
6. Wash wells five times with Wash Buffer making sure each well is filled completely. Remove residual Wash Buffer by gentle tapping or aspiration.
7. Pipette **100 μ L** of **HRP conjugated Detection Antibody** into each well, cover with a plate sealer and incubate **at room temperature (ca.25°C) for 30 minutes**. Discard any unused conjugate.
8. Wash wells five times as same as in step 6.
9. Add **100 μ L** of **Substrate Reagent** to each well and incubate **at room temperature (ca.25°C) for 10-15 minutes**.
- 10 add **100 μ L** of **Stop Solution** to each well in the same order as the previously added Substrate Reagent.
11. Measure absorbance in each well using a spectrophotometric plate reader at dual wavelengths of 450/540 nm. Dual wavelengths of 450/550 or 450/595 nm can also be used. Read the plate at 450



nm if only a single wavelength can be used. Wells must be read within 30 minutes of adding the Stop Solution.

Recommendations

Special considerations when screening activators or inhibitors

In order to estimate the inhibitory effect on CK2 activity in the test chemicals correctly, it is necessary to conduct the control experiment of “Solvent control” at least once for every experiment and “Inhibitor control” at least once for the first experiment, in addition to “Test sample”, as indicated in the following table. When test chemicals cause an inhibitory effect on CK2 activity, the level of A450 is weakened as compared with “Solvent control”. The high level of A450 is not observed in “Inhibitor control” (usually $A450 < 0.3$).

Assay reagents	Test sample	Solvent control	Inhibitor control
Kinase Reaction Buffer	80 μL	80 μL	80 μL
10X Inhibitor or equivalent	10 μL	-	-
Solvent for Inhibitor	-	10 μL	-
10X Heparin (1 μg/mL)*	-	-	10 μL
CK2 Positive Control (2 m unit/μL)* or Your enzyme samples	10 μL	10 μL	10 μL

* See the section “Materials Required but not Provided” above.

1. Following the above table, add the Reagents to each well of the microplate. Finally, initiate reaction by adding 10 μ L of “CK2 Positive Control” or “Your enzyme samples” to each well and mixing thoroughly at room temperature. Cover with plate sealer. Incubate **at 30°C for 30 minutes**.
2. Follow the step 5 to 10 of “Standard Assay” above.

**Special considerations when measuring precise CK2 activity**

In order to measure the activity of CK2 correctly, it is necessary to conduct the control experiment of “Inhibitor control” at least once for every experiment and “ATP minus control” at least once for the first experiment, in addition to “No enzyme control” as indicated in the following table. Although the level of A450 increases in “Test sample” when CK2 enzyme activity is in the sample, the high level of A450 is not observed in “Inhibitor control”, “ATP minus control” and “No enzyme control”.

Assay reagents	Test Sample	Inhibitor control	ATP minus control	Positive control	No enzyme control
Kinase Reaction Buffer	90 µL	80 µL	-	90 µL	90 µL
Kinase Buffer (provided)	-	-	90 µL	-	-
10X Heparin (1 µg/mL)*	-	10 µL	-	-	-
Your enzyme fraction	10 µL	10 µL	10 µL	-	-
CK2 Positive Control (2 m unit/µL)*	-	-	-	10 µL	-
Buffer	-	-	-	-	10 µL

* See the section “Materials Required but not Provided” above.

1. Following the above table, add the Reagents to each well of the microplate. Finally, initiate the reaction by adding 10 µL of “Your enzyme samples” or “CK2 Positive Control” or “Buffer” to each well and mixing thoroughly at room temperature. Cover with plate sealer. Incubate **at 30°C for 30 minutes**.
2. Follow the step 5 to 10 of “Standard Assay” above.



Evaluation of Results

1. Average the absorbance values for the CK2 sample duplicates (positive control) and all experimental sample duplicate values (when applicable). When CK2 Positive Control (20 m units/assay) is included as an internal control for the phosphorylation reaction, the absorbance value should be greater than 1.0 with a background less than 0.15.
2. For screening of purification/chromatography fractions, on graph paper, plot the mean absorbance values for each of the samples on the Y-axis versus the fraction number on the X-axis to determine the location of the eluted, purified CK2.
3. For kinetic analysis, on graph paper, plot the mean absorbance values for each of the time points on the Y-axis versus the time of each reaction (minutes) on the X-axis.

Assay Characteristics

The MBL Research Product **CycLex CK2 Kinase Assay/Inhibitor Screening Kit** has been shown to detect the activity of purified recombinant CK2. The assay shows good linearity of sample response.-

Troubleshooting

1. All samples and standards should be assayed in duplicate, using the protocol described in the **Detailed Protocol**. Incubation times or temperatures significantly different from those specified may give erroneous results.
2. The reaction curve is nearly a straight line if the kinetics of the assay is of the first order. Variations in the protocol can lead to non-linearity of the curve, as can assay kinetics that are other than first order. For a non-linear curve, point to point or quadratic curve fit methods should be used.
3. Poor duplicates, accompanied by elevated values for wells containing no sample, indicate insufficient washing. If all instructions in the **Detailed Protocol** were followed accurately, such results indicate a need for washer maintenance.
4. Overall low signal may indicate that desiccation of the plate has occurred between the final wash and addition of Substrate Reagent. Do not allow the plate to dry out. Add Substrate Reagent immediately after wash.

Reagent Stability

All of the reagents included in the MBL Research Product **CycLex CK2 Kinase Assay/Inhibitor Screening Kit** have been tested for stability. Reagents should not be used beyond the stated expiration date. Upon receipt, kit reagents should be stored at 4°C. Coated assay plates should be stored in the original foil bag sealed by the zip lock and containing a desiccant pack.

Example of Test Results

Fig.1 Dose dependency of recombinant CK2 enzyme reaction




Fig.2 Time course of recombinant CK2 enzyme reaction




Fig.3-1 Effect of TBB (Calbiochem Cat# 218697) on activity of recombinant CK2





Fig.3-2 Effect of Heparin on activity of recombinant CK2





References

1. Lozeman, F.J., Litchfield, D.W., Piening, C., Takio, K., Walsh, K.A. and Krebs, E.G. (1990) Isolation and characterization of human cDNA clones encoding the α and α' subunits of CK2. *Biochemistry* 29, 8436–8447
2. Litchfield, D.W., Lozeman, F.J., Piening, C., Sommercorn, J., Takio, K., Walsh, K.A. and Krebs, E.G. (1990) Subunit structure of CK2 from bovine testis: demonstration that the α and α' subunits are distinct polypeptides. *J. Biol. Chem.* 265, 7638–7644
3. Maridor, G., Park, W., Krek, W. and Nigg, E.A. (1991) CK2. cDNA sequences, developmental expression and tissue distribution of mRNAs for α , α' and β subunits of the chicken enzyme. *J. Biol. Chem.* 266, 2362–2368
4. Munstermann, U., Fritz, G., Seitz, G., Lu, Y.P., Schneider, H.R. and Issinger, O.-G. (1990) CK2 is elevated in solid human tumours and rapidly proliferating non-neoplastic tissue. *Eur. J. Biochem.* 189, 251–257
5. Pepperkok, R., Lorenz, P., Ansorge, W. and Pyerin, W. (1994) CK2 is required for transition of G0/G1, early G1, and G1/S phases of the cell cycle. *J. Biol. Chem.* 269, 6986–6991
6. Landesman-Bollag, E., Romieu-Mourez, R., Song, D.H., Sonenshein, G.E., Cardiff, R.D. and Seldin, D.C. (2001) Protein kinase CK2 in mammary gland tumorigenesis. *Oncogene* 20, 3247–3257
7. Seldin, D.C. and Leder, P. (1995) CK2 α transgene-induce murine lymphoma: relation to theileriosis in cattle. *Science* 267, 894–897
8. Landesman-Bollag, E., Channavajhala, P.L., Cardiff, R.D. and Seldin, D.C. (1998) p53 deficiency and misexpression of protein kinase CK2a collaborate in the development of thymic lymphomas in mice. *Oncogene* 16, 2965–2974
9. Channavajhala, P. and Seldin, D.C. (2002) Functional interaction of protein kinase CK2 and c-Myc in lymphomagenesis. *Oncogene* 21, 5280–5288
10. Sayed, M., Pelech, S., Wong, C., Marotta, A. and Salh, B. (2001) Protein kinase CK2 is involved in G2 arrest and apoptosis following spindle damage in epithelial cells. *Oncogene* 20, 6994–7005
11. Desagher, S., Osen-Sand, A., Montessuit, S., Magnenat, E., Vilbois, F., Hochmann, A., Journot, L., Antonsson, B. and Martinou, J.C. (2001) Phosphorylation of bid by casein kinases I and II regulates its cleavage by caspase 8. *Mol. Cell* 8, 601–611
12. Li, P., Li, J., Muller, E., Otto, A., Dietz, R. and von Harsdorf, R. (2002) Phosphorylation by protein kinase CK2. A signaling switch for the caspase-inhibiting protein ARC. *Mol. Cell* 10, 247–258
13. Wang, H., Davis, A., Yu, S. and Ahmed, K. (2001) Response of cancer cells to molecular interruption of the CK2 signal. *Mol. Cell. Biochem.* 227, 167–174



CK2 Kinase Assay/Inhibitor Screening Kit
User's Manual
For Research Use Only, Not for use in diagnostic procedures



For more information, please visit our website at
<https://ruo.mbl.co.jp/>.

MANUFACTURED BY



URL: <https://ruo.mbl.co.jp>
E-mail: support@mbl.co.jp

CycLex/CircuLex products are supplied for research use only. CycLex/CircuLex products and components thereof may not be resold, modified for resale, or used to manufacture commercial products without prior written approval from MBL. To inquire about licensing for such commercial use, please contact us via email.