

MONOCLONAL ANTIBODY

Anti-CD279 (PD-1) (Human) mAb

Code No.	Clone	Subclass	Quantity	Concentration
D132-3	J110	Mouse IgG1	100 µL	1 mg/mL

BACKGROUND: Human PD-1 (programmed death-1) is a 55 kDa member of the immunoglobulin superfamily that is induced in cells undergoing apoptosis. The PD-1 protein contains an immunoreceptor tyrosine-based inhibitory motif and is expressed predominantly on activated T and B lymphocytes. PD-1 plays a key role in peripheral tolerance and autoimmune disease and is thought to be involved in the maintenance of peripheral self-tolerance by serving as a negative regulator of immune responses. Two novel members of the B7 family have been identified as PD-1 ligands, PD-L1 (B7-H1) and PD-L2 (B7-DC). Evidence reported to date suggests overlapping functions for these two PD-1 ligands and their constitutive expression on some normal tissues and up-regulation on activated antigen-presenting cells.

SOURCE: This antibody was purified from hybridoma (clone J110) supernatant using protein A agarose. This hybridoma was established by fusion of mouse myeloma cell P3U1 with Balb/c mouse splenocyte immunized with the Fc-fusion human PD-1 protein.

FORMULATION: 100 µg IgG in 100 µL volume of PBS containing 50% glycerol, pH 7.2. No preservative is contained.

STORAGE: This antibody solution is stable for one year from the date of purchase when stored at -20°C.

REACTIVITY: This antibody reacts with human PD-1 on Western blotting and Flow cytometry.

APPLICATIONS:

- Western blotting; 1 µg/mL for chemiluminescence detection system
- Immunoprecipitation; Not tested
- Immunohistochemistry; Not tested
- Immunocytochemistry; Not tested
- Flow cytometry; 10-20 µg/mL (final concentration)

Detailed procedure is provided in the following **PROTOCOLS**.

SPECIES CROSS REACTIVITY:

Species	Human	Mouse	Rat
Cell	Transfectant	Not tested	Not tested
Reactivity on WB	+		

INTENDED USE:

For Research Use Only. Not for use in diagnostic procedures.


REFERENCES:

- 1) Bennett, F., *et al.*, *J. Immunol.* **170**, 711-719 (2003)
- 2) Iwai, Y., *et al.*, *Immunol. Lett.* **83**, 215-220 (2002)

Clone J110 is used in these references.

RELATED PRODUCTS:

- D132-3 Anti-CD279 (PD-1) (Human) mAb (J110)
- D132-4 Anti-CD279 (PD-1) (Human) mAb-FITC (J110)
- D132-5 Anti-CD279 (PD-1) (Human) mAb-PE (J110)
- D133-3 Anti-CD279 (PD-1) (Human) mAb (J105)
- D133-5 Anti-CD279 (PD-1) (Human) mAb-PE (J105)
- D092-3 Anti-CD274 (PD-L1) (Human) mAb (MIH3)
- D092-6 Anti-CD274 (PD-L1) (Human) mAb-Biotin (MIH3)
- D230-3 Anti-CD274 (PD-L1) (Human) mAb (27A2)
- D231-3 Anti-CD273 (PD-L2) (Mouse) mAb (54-1)



Western blot analysis of PD-1 expression in X63 cells (1) and human PD-1 transfected X63 cells (2) using D132-3.

PROTOCOLS:

SDS-PAGE & Western Blotting

- 1) Wash the cells 3 times with PBS and suspend with 10 volume of cold Lysis buffer (50 mM Tris-HCl, pH 7.2, 250 mM NaCl, 0.1% NP-40, 2 mM EDTA, 10% glycerol) containing appropriate protease inhibitors. Incubate it at 4°C with rotating for 30 minutes, then sonicate briefly (up to 10 seconds).
- 2) Centrifuge the tube at 12,000 x g for 10 minutes at 4°C and transfer the supernatant to another tube. Measure the

protein concentration of the supernatant and add the cold Lysis buffer to make 8 mg/mL solution.

- 3) Mix the sample with equal volume of Laemmli's sample buffer.
- 4) Boil the samples for 3 minutes and centrifuge. Load 10 μ L of the sample per lane in a 1 mm thick SDS-polyacrylamide gel for electrophoresis.
- 5) Blot the protein to a polyvinylidene difluoride (PVDF) membrane at 1 mA/cm² for 1 hour in a semi-dry transfer system (Transfer Buffer: 25 mM Tris, 190 mM glycine, 20% MeOH). See the manufacture's manual for precise transfer procedure.
- 6) To reduce nonspecific binding, soak the membrane in 10% skimmed milk (in PBS, pH 7.2) for 1 hour at room temperature, or overnight at 4°C.
- 7) Incubate the membrane with primary antibody diluted with PBS, pH 7.2 containing 1% skimmed milk as suggest in the **APPLICATIONS** for 1 hour at room temperature. (The concentration of antibody will depend on condition.)
- 8) Wash the membrane with PBS-T [0.05% Tween-20 in PBS] (5 minutes x 3 times).
- 9) Incubate the membrane with the 1:10,000 Anti-IgG (Mouse) pAb-HRP (MBL; code no. 330) diluted with 1% skimmed milk (in PBS, pH 7.2) for 1 hour at room temperature.
- 10) Wash the membrane with PBS-T (10 minutes x 3 times).
- 11) Wipe excess buffer on the membrane, then incubate it with appropriate chemiluminescence reagent for 1 minute.
- 12) Remove extra reagent from the membrane by dabbing with paper towel, and seal it in plastic wrap.
- 13) Expose to an X-ray film in a dark room for 3 minutes.
- 14) Develop the film as usual. The condition for exposure and development may vary.

(Positive control for Western blotting; Transfectant)

Flow cytometric analysis for floating cells


We usually use Fisher tubes or equivalents as reaction tubes for all step described below.

- 1) Wash the cells 3 times with washing buffer [PBS containing 2% fetal calf serum (FCS) and 0.1% NaN₃].
- 2) Resuspend the cells with washing buffer (5x10⁶ cells/mL).
- 3) Add 50 μ L of the cell suspension into each tube, and centrifuge at 500 x g for 1 minute at room temperature (20~25°C). Remove supernatant by careful aspiration.
- 4) Add 10 μ L of normal goat serum containing 1 mg/mL normal human IgG and 0.1% NaN₃ to the cell pellet after tapping. Mix well and incubate for 5 minutes at room temperature.
- 5) Add 40 μ L of the primary antibody at the concentration of as suggested in the **APPLICATIONS** diluted with the washing buffer. Mix well and incubate for 30 minutes at room temperature.
- 6) Add 1 mL of the washing buffer followed by centrifugation at 500 x g for 1 minute at room


temperature. Remove supernatant by careful aspiration.

- 7) Add 30 μ L of 1:100 Anti-IgG (Mouse) pAb-FITC (MBL; code no. IM-0819) diluted with the washing buffer. Mix well and incubate for 15 minutes at room temperature.
- 8) Add 1 mL of the washing buffer followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 9) Resuspend the cells with 500 μ L of the washing buffer and analyze by a flow cytometer.

(Positive controls for Flow cytometry; Transfectant, peripheral blood lymphocyte)



Flow cytometric analysis of PD-1 expression on human PD-1 transfected X63 cells (left) and the parental cells (right). Open histogram indicates the reaction of isotypic control to the cells. Shaded histograms indicate the reaction of D132-3 to the cells.



Flow cytometric analysis of PD-1 expression on human peripheral blood lymphocytes (1), peripheral blood monocytes (2) and peripheral blood neutrophils (3). Dot line histograms indicate the reaction of isotypic control to the cells. Solid line histograms indicate the reaction of D132-3 to the cells.