CODE No. 7630

Revised: December, 2021 ver. 7

For Research Use Only. Not for use in diagnostic procedures.

Quantitative test kit for Human EDN

EDN ELISA Kit

CODE No. 7630 96 wells

CONTENTS

1. Intended Use -----	1
2. Summary and Explanation-----	1
3. Principle -----	1
4. Species cross reactivity-----	2
5. Materials provided -----	2
6. Storage and Stability -----	2
7. Materials and equipment required -----	3
8. Precautions -----	3
9. Procedure -----	4
10. Performance Characteristics -----	8
11. References -----	10

CODE No. 7630

Before use, thoroughly read all Instructions.

Intended Use

The **EDN ELISA Kit** is based on sandwich ELISA and is capable of measuring Human EDN.

For Research Use Only. Not for use in diagnostic procedures.

Summary and Explanation

Upon activation, eosinophils secrete several cationic proteins including major basic protein (MBP), eosinophil cationic protein (ECP), eosinophil peroxidase (EPO) and eosinophil-derived neurotoxin (EDN). EDN, also called eosinophil protein X (EPX), is a glycosylated single-chain protein with a molecular weight 18-21 kDa and is a member of the ribonuclease superfamily along with ECP¹¹⁻¹³). Although EDN has high sequence homology to ECP, EDN has 100-fold greater ribonuclease activity than ECP and neurotoxicity but not cytotoxicity^{9), 10)}. Eosinophil activation has been found to play an important role in inflammatory processes in allergic diseases. EDN can be a marker of eosinophil activation and degranulation.

The “EDN ELISA Kit” is a useful reagent for specifically measuring Human EDN with high sensitivity by ELISA.

Principle

The EDN ELISA Kit measures Human EDN by sandwich ELISA. This ELISA detects Human EDN with a minimum detection limit of 0.62 ng/mL and does not cross-react with ECP.

Samples to be measured or standards are incubated in the wells coated with anti-human EDN monoclonal antibody. After washing, the peroxidase conjugated anti-human EDN polyclonal antibody is added into the microwells and incubated. After another washing, the peroxidase substrate is mixed with the chromogen and allowed to incubate for an additional period of time. An acid solution is added to each microwell to terminate the enzyme reaction and to stabilize the color development. The optical density (O.D.) of each microwell is measured at 450 nm using a microplate reader. The concentration of Human EDN is calibrated from a dose response curve based on the reference standards.

Species cross reactivity

Species	Human	Mouse	Rat
Samples	serum, plasma, urine	serum	Not Tested
Reactivity	+	-	

Materials provided

Each kit contains;

Name	Materials	Quantity (1 kit)
Microwells	microwell strips coated with anti-human EDN antibody	8-well strip × 12 strips
EDN standard	The standards of 7 types of concentrations (ready-for-use)	1 mL × 7 vials
Assay diluent	Buffer for diluting samples* (ready-for-use)	50 mL × 1 bottle
Conjugate	Peroxidase conjugated anti-human EDN antibody (ready-for-use)	12 mL × 1 bottle
Wash concentrate	Buffer for washing microwells (×10)	100 mL × 1 bottle
Substrate	TMB/H ₂ O ₂ solution (ready-for-use)	20 mL × 1 bottle
Stop solution	0.5 M H ₂ SO ₄ (ready-for-use) (irritant)	20 mL × 1 bottle
Positive control		0.2 mL × 1 vial
Negative control		0.2 mL × 1 vial

* Use the **Assay diluent** as 0 ng/mL control for the EDN standard. During storage, some precipitation or turbidity may appear. It does not affect the assay results. Dissolve the Assay diluent completely using a vortex mixer.

Storage and Stability

All kit components must be stored at 2-8°C. All reagents are stable for 18 months after manufacturing when stored at the indicated conditions.

Materials and equipment required

- Microplate reader (450 nm)
- Automatic plate washer or wash bottle
- Adjustable micropipette
- Multichannel micropipette
- 96-well polyvinyl plate
- Microplate holder
- Uncoated microwell strips (for using with an automatic plate washer)
- Disposable reagent vessels
- Paper towels
- Distilled water

Precautions

1. Allow all the components to come to room temperature (20-25°C) before use.
2. **Microwells** that are not used immediately should be returned to the ziplock pouch, which must be carefully resealed to avoid moisture absorption.
3. Fresh samples should be used. Aliquot each sample and store below -20°C if necessary. Avoid repeated freezing and thawing. Never store the samples at 4°C, as samples are adversely affected by storage at this temperature.
4. When **Assay diluent** is stored at 2-8°C, some precipitation or turbidity may appear. However it does not affect the reagent efficiency. Dissolve the **Assay diluent** completely before use.
5. **Assay diluent** and **EDN Standard** contain sodium azide (0.09%) as a preservative. Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush drain with plenty of water when disposing materials containing azide.
6. **Stop Solution** is 0.5 M sulfuric acid. As it is corrosive product, protect eyes and skin and handle with care.
7. This kit is intended for research use only. Not for use in diagnostic procedures.

Procedure

◆ Preparation of Reagents

1. Wash solution

Prepare “**Wash solution**” by diluting 100 mL of **Wash concentrate** with 900 mL of distilled water. This **Wash solution** is stable for 2 weeks at 4°C.

2. Other reagents

EDN standard, Conjugate, Assay diluent, Substrate and **Stop solution** are ready for use.

◆ Preparation of samples

1. Sampling

Serum levels of EDN can vary significantly due to *in vitro* effects after sampling. These effects may be due to EDN released during coagulation. Therefore, each laboratory is recommended to establish its own sampling procedure. The following sampling procedures have been reported.

1) Serum

Blood for measurement of EDN was collected using Venoject®-II Tube AUTOSEP tube (TERUMO, Japan; cat# VP-AS109K). The blood was allowed to clot at room temperature (20-25°C) for exactly 60 ± 10 minutes. Thereafter, the serum was separated by centrifugation twice for 10 minutes at $1,350 \times g$ at 4°C, and then transferred to a new plastic test tube.

2) Plasma

Blood for measurement of EDN was collected using Venoject®-II Tube (TERUMO, Japan; cat# VP-DK050K). Within 30 minutes of blood collection, the plasma was separated by centrifugation twice for 10 minutes at $1,350 \times g$ at 4°C, and then transferred to a new plastic test tube.

3) Urine

Within 30 minutes of urine collection, the urine was separated by centrifugation, twice for 10 minutes at $1,350 \times g$ at 4°C, and then transferred to a new plastic test tube.

- 24-hour total urine sample; Daily urinary excretion of EDN was presented as mg/day.
- When the 24-hour total urine was not available, urine was collected from each individual as a spot sample. Urinary creatinine was also quantitated and urinary EDN results were presented as $\mu\text{g}/\text{mmol}$ creatinine.

2. Dilution

1) Dilute samples

Dilute each sample with **Assay diluent**.

Example 1: human serum or plasma

Dilute 1:5 with **Assay diluent** by adding 50 μL of sample to 200 μL of **Assay diluent**.

Example 2: human urine

Dilute 1:51 with **Assay diluent** by adding 10 μL of sample to 500 μL of **Assay diluent**.

2) Dilute **Positive control** and **Negative control**

Dilute each control 1:5 with **Assay diluent** by Adding 50 μL of each control to 200 μL of **Assay diluent**.

3. Storage

Fresh samples should be used. Aliquot each sample into a new plastic tube, and store below -20°C if necessary. Avoid repeated freezing and thawing.

◆ Assay procedure

Perform all assays in duplicate.

STEP 1: Sample incubation

- 1) Add 150 μL of prepared samples, **EDN standard** (ready-for-use), and controls to a transient 96-well polyvinyl plate as the same order of actual assay run. Then, transfer 100 μL of each sample to **Microwells** simultaneously using a multichannel micropipette.

The reaction starts immediately upon pipetting to **Microwells**. This transfer should be completed as quickly as possible.

- 2) Incubate for 60 minutes at room temperature ($20-25^{\circ}\text{C}$).

STEP 2: Washing

Manual wash: Aspirate or discard the well contents, then fill each well with **Wash solution** using a wash bottle. Aspirate or discard the **Wash solution** in the wells. Repeat this step another 3 times. Tap the plate on a paper towel to remove any remaining **Wash solution**.

Autowasher: Wash 4 times with **Wash solution**. Tap the plate on a paper towel to remove any remaining **Wash solution**.

* Each laboratory is recommended to confirm its own appropriate washing times and set-up.

* **Wash solution** should be used at room temperature ($20-25^{\circ}\text{C}$).

STEP 3: Incubation with Conjugate

- 1) Pour **Conjugate** (ready-for-use) into the vessel. After completely removing any remaining wash solution, pipette 100 μL of **Conjugate** into each well with a multichannel micropipette.
- 2) Incubate for 60 minutes at room temperature (20-25°C).

STEP 4: Washing

Wash the wells again following the **STEP 2** procedure.

STEP 5: Substrate incubation

- 1) After the washing, pour **Substrate** (ready-for-use) into vessel. Add 100 μL of **Substrate** to each well.

The reaction starts immediately upon pipetting to the well. This transfer should be completed as quickly as possible by using multichannel micropipette.

- * **Substrate** should be used at room temperature (20-25°C).
- * Use disposable new pipettes and vessel, as **Substrate** is easily oxidized by metal ions and contaminated by microbes.
- * Once **Substrate** is poured into the vessel, do not return to the bottle.

- 2) Incubate for 10 minutes at room temperature (20-25°C).

STEP 6: Stopping reaction

Pour **Stop solution** (ready-for-use) into the vessel. Add 100 μL of **Stop solution** to each well with a multichannel micropipette.

◆ Reading

Read the absorbance of each well at 450 nm using a plate reader. If a dual wavelength plate reader is available, set the test wavelength at 450 nm and the reference at 620 nm.


- * Reading should be done within 30 minutes after stopping the reaction.
- * Ensure that the back of the plate is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before reading.

◆ **Calculation of results**

- 1) Calculate the mean absorbance value of each **EDN standard**.
- 2) Plot on the semi-log graph paper and construct a calibration curve (absorbance on the vertical axis, concentration (ng/mL) on the horizontal axis).
- 3) Read the ng/mL corresponding to the absorbance of sample from the standard curve.
Report the EDN concentration of samples by multiplying dilution factor. (e.g. human serum: ×5)

* If absorbance of sample exceeds the absorbance reading for the highest **EDN standard**, dilute the sample and measure again.

◆ **Example of calibration curve**

◆ **Concentration of Human EDN in normal human sample**

Serum, plasma and urine samples from healthy donors were measured with the EDN ELISA Kit. The results were as follows. (Calculated as described in ◆ **Calculation of results**.)

	serum	plasma	urine
Sample (n)	52	52	50
Maximum	66.4 ng/mL	49.8 ng/mL	164.2 µg EDN/mmol creatinine
Minimum	8.3 ng/mL	6.2 ng/mL	26.7 µg EDN/mmol creatinine
Mean	26.4 ng/mL	18.1 ng/mL	81.8 µg EDN/mmol creatinine
SD	13.3	9.6	36.2
Mean +3 SD	66.3 ng/mL	46.9 ng/mL	190.4 µg EDN/mmol creatinine

Performance Characteristics

◆ Sensitivity

The sensitivity of the assay is 0.62 ng/mL.

◆ Reproducibility

1. Intra-assay

Intra-assay reproducibility was determined by assaying the sera 9 times.

EDN concentrations of the serum samples were calculated as described in

◆ Calculation of results.

Sample	serum 1	serum 2	serum 3
Number of determinations	9	9	9
Mean (ng/mL)	1.9	3.5	16.0
C.V. (%)	3.6	2.8	2.6

2. Inter-assay

Inter-assay reproducibility was determined by 4 independent “Intra-assay”.

EDN concentrations of the serum samples were calculated as described in

◆ Calculation of results.

Sample	serum 1	serum 2	serum 3
Number of determinations	4	4	4
Mean (ng/mL)	3.8	4.2	10.3
C.V. (%)	6.5	7.2	9.4

* From 4 replicates of each serum sample in 4 separate assays.

◆ **Recovery test**

Purified Human EDN was added to samples at different concentrations.

EDN concentrations of the serum samples were calculated as described in

◆ **Calculation of results.**

serum 1

(A) additional EDN (ng/mL)	EDN concentration observed (ng/mL)	(B) recovery (ng/mL)	(B/A) recovery (%)
0.0	4.5	-	-
2.1	6.6	2.1	100
4.1	8.6	4.1	100
8.2	12.6	8.1	99

serum 2

(A) additional EDN (ng/mL)	EDN concentration observed (ng/mL)	(B) recovery (ng/mL)	(B/A) recovery (%)
0.0	5.6	-	-
2.1	7.5	1.9	90
4.1	9.1	3.5	85
8.2	14.6	9.0	110


serum 3

(A) additional EDN (ng/mL)	EDN concentration observed (ng/mL)	(B) recovery (ng/mL)	(B/A) recovery (%)
0.0	14.3	-	-
2.1	16.2	1.9	90
4.1	18.1	3.8	93
8.2	21.7	7.4	90

◆ Dilution test

Serum samples were diluted with **Assay diluent**.

Each serum was serially diluted from 1/1 to 1/128 with **Assay diluent**, then EDN concentrations of the serum samples were calculated as described in ◆ **Calculation of results**.

References

- 1) Månsson, A., *et al.*, *J. Leu. Biol.* **85**, 719-727 (2009)
- 2) Gaudreault, É., *et al.*, *J. Immunol.* **180**, 6211-6221 (2008)
- 3) Mathur, S. K., *et al.*, *CHEST* **133**, 412-419 (2008)
- 4) Munitz, A., *et al.*, *J. Immunol.* **177**, 77-83 (2006)
- 5) Morokata, T., *et al.*, *J. Pharmacology and Experimental Therapeutics* **317**, 244-250 (2006)
- 6) Decot, V., *et al.*, *J. Immunol.* **174**, 628-635 (2005)
- 7) Thomas, L. L., *et al.*, *J. Immunol.* **169**, 993-999 (2002)
- 8) Nagase, H., *et al.*, *J. Immunol.* **164**, 5935-5943 (2000)
- 9) Slifman, N. R., *et al.*, *J. Immunol.* **143**, 2317-2322 (1989)
- 10) Gullberg, U., *et al.*, *Biochem. Biophys. Res. Commun.* **139**, 1239-1242 (1986)
- 11) Slifman, N. R., *et al.*, *J. Immunol.* **137**, 2913-2917 (1986)
- 12) Peterson, C. G. B., *et al.*, *Immunology* **50**, 19-26 (1983)
- 13) Durack, D.T., *et al.*, *PNAS.* **78**, 5165-5169 (1981)

CODE No. 7630

As this product has been used in many researches, these references are a part of such reports.
Please visit our website at <https://ruo.mbl.co.jp/>.

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