# Anti-Integrin $\alpha v \beta 6$ ELISA Kit

CODE 5288

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**For Research Use Only, Not for use in diagnostic procedures**

Enzyme immunoassay for the detection of anti-Integrin  $\alpha v \beta 6$  antibody in serum

# Anti-Integrin $\alpha v \beta 6$ ELISA Kit

Intended Use.....	1
Storage.....	1
Introduction.....	2
Principle of the Assay.....	2
Materials Provided.....	2
Detailed Protocol.....	3
Warnings and Precautions.....	5
Assay Characteristics.....	7
References.....	9

## Intended Use

The MBL Research Product **Anti-Integrin  $\alpha v \beta 6$  ELISA Kit** is intended to be used for quantitative measurement of anti-integrin  $\alpha v \beta 6$  antibody in serum samples.

The Standard material and Reaction buffer contains 0.09% sodium azide as a preservative. If these reagents get into the eyes or mouth, or on the skin, first-aid measures should be taken, such as flushing with water, in addition to seeking medical attention if deemed necessary.

Carefully handle the Stop solution and avoid any contact with the eyes, mouth, or skin because this reagent contains sulfuric acid. If contact occurs, immediately flush the affected area with copious amounts of water and seek medical attention if deemed necessary.

**This assay kit is for research use only and is not intended for diagnostic or therapeutic procedures.**

## Storage

- Upon receipt, store all components at 2–8°C.
- The expiration date of the reagents is indicated on the outer packaging.

---

## Introduction

---

Integrins are heterodimeric cell surface receptors involved in cell-extracellular matrix and cell-cell adhesion. In addition, integrins convey two-way signaling between cells and their surroundings and participate in the regulation of a variety of cellular processes, including cell proliferation, migration, survival, tissue invasion, and innate immunity. Mammalian cells have 18 distinct integrin  $\alpha$  subunits and 8  $\beta$  subunits, which can combine in various ways to form 24 different receptor heterodimers. Integrin  $\alpha\beta6$  is expressed exclusively by epithelial cells and functions as a receptor for extracellular matrix proteins such as fibronectin. Furthermore, integrin  $\alpha\beta6$  expressed in the intestinal epithelium is involved in suppression of inflammation and prevention of pathogen and parasite infection [1]. Anti-integrin  $\alpha\beta6$  antibodies have been detected in the vast majority of the patients with ulcerative colitis, and its titers have been reported to correlate with disease activity [2].

---

## Principle of the Assay

---

The Anti-Integrin  $\alpha\beta6$  ELISA Kit readily allows measurement of anti-integrin  $\alpha\beta6$  antibody levels in human serum. Calibration controls and diluted serum samples are added to ELISA plates coated with integrin  $\alpha\beta6$ . After washing away unbound serum proteins, anti-human IgG conjugated with horseradish peroxidase is added and the plate is incubated (Conjugate incubation). After a wash step, the chromogenic substrate for peroxidase is added and the plate is incubated (Substrate incubation). Finally, an acidic solution is added to stop the enzyme reaction and consequent color development. Anti-integrin  $\alpha\beta6$  IgG concentrations can then be determined by comparing the absorbance at 450 nm of samples relative to those of calibration controls.

---

## Materials Provided

---

Components	Package	Quantity
$\alpha\beta6$ coated microplate	8-well x 12 strips	1 pack
Standard material	50 $\mu$ L	1 vial
Reaction buffer	50 mL	1 bottle
Wash concentrate (10x)	100 mL	1 bottle
Conjugate diluent	20 mL	1 bottle
HRP conjugated antibody	50 $\mu$ L	1 vial
Substrate reagent	20 mL	1 bottle
Stop solution	20 mL	1 bottle
Plate seal	-	3 sheets

- Reagents from different lot numbers should not be combined.
- All components should be stored at 2–8°C.
- A plate seal is provided to prevent dust and evaporation of the reaction solution during the first and second reactions.

## Detailed Protocol

### Preparation of required items

- Microplate reader
- Plate washer or wash bottle
- Micropipette
- Multichannel micropipette
- Reagent reservoir
- ddH<sub>2</sub>O
- 15 mL centrifuge tube

### Preparation of Working Solutions

Allow all the reagents to reach **room temperature (approx. 20–25°C)** before use.

1. Prepare the wash solution by adding **100 mL** of the **Wash concentrate (10x)** to **900 mL** of deionized (distilled) water (**ddH<sub>2</sub>O**). The wash solution should be at approximately **25°C** when used. The 1X-diluted wash solution is stable for 2 weeks at 2–8°C.
2. Dilute the **HRP conjugated antibody** with the **Conjugate diluent**. Dispense **12 mL** of the **Conjugate diluent** into a 15 mL centrifuge tube. Add **30 µL** of **HRP conjugated antibody** and mix well. Use the prepared enzyme-labeled antibody solution at **room temperature (approx. 20–25°C)**.

### Assay Procedure

#### (1) Standard solution and sample preparation

Use the **Standard material** to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 200.0 U/mL standard (Cal. 7) serves as the highest standard. The **Reaction buffer** serves as the zero standard (Blank).

	Volume of Standard	Assay Diluent (red)	Concentration (U/mL)
Cal. 7	5 µL of 101X Standard solution	500 µL	200.0
Cal. 6	250 µL of Cal. 7	250 µL	100.0
Cal. 5	250 µL of Cal. 6	250 µL	50.0
Cal. 4	250 µL of Cal. 5	250 µL	25.0
Cal. 3	250 µL of Cal. 4	250 µL	12.5
Cal. 2	250 µL of Cal. 3	250 µL	6.3
Cal. 1	250 µL of Cal. 2	250 µL	3.1
Blank	—	250 µL	0.0

Dilute each serum sample 1:101 by adding 5 µL of serum to 500 µL of **Reaction buffer**.

#### (2) 1<sup>st</sup> Reaction

- 1) Remove the appropriate number of microtiter wells from the  **$\alpha\text{v}\beta\text{6}$  coated microplate** and place them into the well holder. Return any unused wells to the foil pouch, refold, seal with tape, and store at 4°C.
- 3) Pipette **100 µL** of **Calibrators (Blank and Cal. 1–7)** and **diluted samples** into the appropriate wells.

Note: The reaction starts as soon as the sample is added to the  **$\alpha\text{v}\beta\text{6}$  coated microplate**, so the pipetting into the plate needs to be performed quickly and without interruptions at every step of the test procedure.

- 4) Cover the wells with a plate seal and incubate the plate for 60 minutes at **room temperature**

**(approx. 20–25°C).**

(3) Washing

- 1) Remove the solution in the wells.
- 2) Wash four times by filling each well with 1X wash solution (> 350  $\mu\text{L}$ ) using a squirt bottle, multichannel micropipette, manifold dispenser, or microplate washer.

Note-1: The wash solution should be used at room temperature (approx. 20–25°C).

Note-2: When using an automatic washer, the appropriate washing conditions should be established by the individual user.

(4) 2<sup>nd</sup> Reaction

- 1) Remove all of the wash buffer from the wells.
- 2) Add 100  $\mu\text{L}$  of the **1X HRP conjugated antibody** to each well using a multichannel micropipette.

Note: The **HRP conjugated antibody** should not be returned to the reagent bottle in order to avoid contamination between reagents.

- 3) Cover the wells with the plate seal and incubate the plate for 60 minutes at **room temperature (approx. 20–25°C).**

(5) Washing

- 1) Remove the solution in the wells.
- 2) Wash four times, as in step (3)-2).

(6) Enzyme reaction

- 1) Fully remove the wash buffer in the wells.
- 2) Add **100  $\mu\text{L}$  of Substrate reagent** to each well using a multichannel micropipette.

Note-1: Avoid exposing the microtiter plate to direct sunlight.

Note-2: Use a fresh reservoir for the Substrate reagent to avoid contamination of the reagents.

Note-3: The Substrate reagent should not be returned to the bottle after using it.

- 3) Incubate for 20 minutes at **room temperature (approx. 20–25°C).**

(7) Stop Reaction

Add **100  $\mu\text{L}$  of Stop solution** to each well using a multichannel micropipette in the same order as in step (6).

(8) Reading

Measure the absorbance in each well using a spectrophotometric microplate reader at dual wavelengths of 450 and 620 nm. Read the microplate at 450 nm if only a single wavelength can be used. Wells must be read within 30 minutes of adding the Stop solution.

Note: Complete removal of the liquid at each step is essential for good performance. After the last wash, remove any remaining wash solution by aspiration or decanting. Invert and tap the plate on clean paper towels to remove any residual traces of fluid.

(9) Calculation

Generate a calibration curve using the OD value of Cal. 1–7 and Blank at 450 nm (y-axis) and the antibody concentration (x-axis). The calibration curve fits best to a sigmoidal four-parameter logistic equation. If an OD value of a sample is higher than that of the highest concentration calibrator (Cal. 7), re-testing with further diluted samples is recommended.

Each assay result should meet the following criteria.

OD value (450nm) of Blank:  $\leq 0.100$

OD value (450nm) of Cal. 7:  $\geq 1.0000$

---

## **Warnings and Precautions**

---

1. All patient samples should be treated as if they pose a risk of infection, such as HIV, HBV, HCV, etc.
2. If these reagents get into the eyes or mouth, or on the skin, first-aid measures should be taken, such as flushing with water, in addition to seeking medical attention if deemed necessary.
3. Carefully handle the Stop solution and avoid any contact with the eyes, mouth, or skin because this reagent contains sulfuric acid. If contact occurs, immediately flush the affected area with copious amounts of water and seek medical attention if deemed necessary.
4. This assay kit is for research use only and is not intended for use in diagnostic or therapeutic procedures.
5. Do not use any of the kit components beyond the stated expiration dates.
6. Do not mix reagents from different lot numbers.
7. All reagents must be brought to room temperature (approx. 20–25°C) before starting the assay.
8. Do not expose the kit components to direct sunlight during storage or when in use during the assay.
9. Avoid microbial contamination of all reagents and samples.
10. Strict adherence to the specified assay instructions (incubation temperature, incubation time, and sample dilution) is essential.
11. The wells must be properly rinsed with 1X wash solution to avoid false-positive or false-negative results.
12. Do not allow the microplate to dry out during the assay.
13. Avoid cross-contamination of microtiter wells. Mix reagents gently so as not to form bubbles in the liquid or allow droplets to splash out.
14. Use disposable pipettes when transferring the Substrate reagent from the reagent bottle to the reservoir because the Substrate is readily oxidized by metal ions.  
Do not return the Substrate transferred to the reservoir back to the reagent bottle.
15. All breakaway microtiter wells that are not immediately required should be returned to the pouch with desiccant and must be carefully resealed to avoid moisture from being absorbed.
16. The Reaction buffer may form a precipitate. Before use in the test, bring the buffer to room temperature (approx. 20–25°C) and thoroughly mix it to dissolve any precipitate. This will not affect the test results.
17. The Wash concentrate (10x) may look turbid owing to crystal formation at 2-8°C, but this does not cause inconsistent results. Prior to preparation of the 1X wash buffer working solution, bring it to

**For Research Use Only, Not for use in diagnostic procedures**

25°C and mix it thoroughly to dissolve any material that has formed a precipitate.

18. In accordance with the principles of good laboratory practice, it is strongly recommended that all clinical specimens and assay materials, including the antigen microtiter wells, should be considered potentially infectious and handled and disposed of with all necessary precautions. Disposable implements (pipettes, tips, tubes, etc.) used in the test should be soaked in a final concentration of 2% glutaraldehyde solution or 0.1% sodium hypochlorite solution for at least one hour or autoclaved at 121°C for 20 minutes prior to being discarded.

19. If the serum samples etc. splatter, wear gloves and wipe the contaminated surfaces with a paper towel and sodium hypochlorite.

20. Some of the reagents contain sodium azide as a preservative. Azide reacts with copper or lead in the plumbing system to form explosive metal azides. Therefore, always flush with ample amounts of water when draining materials containing azide. Sodium azide is toxic if ingested.

## Assay Characteristics

### Accuracy

The ratio of the measured value to the expected value for 9 samples with known concentrations were within 0.80–1.20.

	Expected value (U/mL)	1st		2nd		3rd	
		Measured value (U/mL)	Measured/Expected	Measured value (U/mL)	Measured/Expected	Measured value (U/mL)	Measured/Expected
<b>Sample 1</b>	115.4	125.5	1.09	118.7	1.03	114.5	0.99
<b>Sample 2</b>	89.8	99.2	1.11	91.6	1.02	90.8	1.01
<b>Sample 3</b>	75.4	79.0	1.05	74.8	0.99	74.6	0.99
<b>Sample 4</b>	42.6	47.5	1.11	44.0	1.03	42.9	1.01
<b>Sample 5</b>	35.4	38.2	1.08	36.4	1.03	36.8	1.04
<b>Sample 6</b>	22.6	25.2	1.11	23.7	1.05	20.5	0.91
<b>Sample 7</b>	4.9	5.1	1.03	5.0	1.01	5.3	1.07
<b>Sample 8</b>	4.6	4.9	1.08	4.4	0.96	4.6	1.01
<b>Sample 9</b>	4.4	4.7	1.07	4.2	0.95	4.8	1.09

### Reproducibility


Intra-assay reproducibility was determined by testing 9 samples 8 times. The coefficient of variation (CV) (%) values for the reproducibility and repeatability were below 10% for each sample.


	Sample No.								
	1	2	3	4	5	6	7	8	9
<b>1</b>	134.4	104.1	86.6	48.3	40.0	26.5	5.3	5.3	5.4
<b>2</b>	123.1	101.6	83.9	47.9	40.3	25.3	5.2	5.1	4.9
<b>3</b>	126.9	98.1	78.5	48.0	38.1	25.5	5.0	4.5	4.5
<b>4</b>	123.4	96.8	79.4	45.7	38.3	25.0	5.0	4.8	4.6
<b>5</b>	123.5	95.1	77.9	48.0	38.1	24.5	5.0	4.8	4.5
<b>6</b>	122.4	100.2	78.1	47.6	35.9	24.3	5.1	4.9	4.6
<b>7</b>	121.2	96.5	78.6	46.9	36.5	24.6	5.0	4.9	4.6
<b>8 (U/mL)</b>	129.5	101.9	70.1	47.3	38.8	25.5	5.1	5.1	4.7
<b>mean (U/mL)</b>	125.5	99.3	79.1	47.5	38.2	25.2	5.1	4.9	4.7
<b>SD</b>	4.5	3.1	4.8	0.8	1.5	0.7	0.1	0.2	0.3
<b>CV (%)</b>	3.6	3.2	6.1	1.7	4.0	2.9	2.3	4.8	6.1

### Sensitivity

The limit of quantitation (LOQ) of anti-integrin  $\alpha\beta6$  antibodies was 0.4 U/mL.

A regression equation was developed from the mean value of the measurements and the CV of the dilution series of a low level of anti-integrin  $\alpha\beta6$  antibodies in octuplicate. The LOQ was calculated as the value at CV of 10%.

**Typical Calibration Curve**

---

**References**

---

1. Koivisto L, Bi J, Häkkinen L, Larjava H. Integrin  $\alpha\beta6$ : Structure, function and role in health and disease. *Int J Biochem Cell Biol.* 2018;99:186-96.
2. Kuwada T, Shiokawa M, Kodama Y, Ota S, Kakiuchi N, Nannya Y, et al. Identification of an Anti-Integrin  $\alpha\beta6$  Autoantibody in Patients With Ulcerative Colitis. *Gastroenterology.* 2021;160(7):2383-94.e21.

MANUFACTURED BY

**MBL** A JSR Life Sciences Company  
MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.  
URL: <https://ivd.mbl.co.jp/>  
E-mail: [kensa@mbl.co.jp](mailto:kensa@mbl.co.jp)